

dr. aktive[®]

ACHILLES RUPTURE NIGHT SPLINT

What it treats:

- ✓ Specifically designed for Achilles tendon rupture (complete tear of the tendon).
- ✓ The Achilles Rupture Night Splint ensures the foot is kept in plantar flexion (foot pointed down) for the first 3 months after an Achilles tendon tear.

When to use:

- ✓ Use instead of your orthopaedic boot only when resting (sitting or lying down).



Reinvented recovery:

Currently patients wear a large, outdoor, orthopaedic boot day and night for their 10 week rehabilitation. Patients find these boots hot, uncomfortable, dirty and heavy. The Achilles Rupture Night Splint provides a solution to these issues to support patient comfort and recovery.

"Setting a new standard in Achilles Tendon rehabilitation"

Foot and Ankle Surgeon



Features & Benefits:

- Sleep and rest comfortably
- Recovery becomes more bearable
- Cool and lightweight for the patient
- Contoured for the patient's foot
- Better quality of sleep
- Clean indoor splint



dr. aktive[®]

ACHILLES RUPTURE NIGHT SPLINT

ACHILLES TENDON NIGHT ORTHOSIS			SHOE SIZE	
LEFT	RIGHT	SIZE	MENS	WOMENS
DAATRNS-SL	DAATRNS-SR	Small	<8.5	<8.5
DAATRNS-LL	DAATRNS-LR	Large	8.5 and up	8.5 and up

Contact us now to place your order or for more information:

www.kinetecuk.com

enquiries@kinetecuk.com

+44 (0)3330 347400

