


T55 - Montage complet sur lit Stryker SV1 - SV2


T55 - Complete set on Stryker SV1 - SV2 bed

Adaptation spéciale pour lit :
Special Adaptation for bed:

stryker®
SV1 - SV2



Dimensions en mm
(*Lit en position basse)
Sizes in mm
(*Bed in low position)



4655005188		T55 - Montage complet sur lit Stryker SV1 - SV2	T 55 - Complete set on Stryker SV1 - SV2 bed
	Q.	Le cadre est composé de :	The frame consists of these parts:
4655005162	1	PRTT55 - Portique tête de lit Stryker SV1 - SV2	PRTT55 - Head porch bar for Stryker SV1 - SV2
4655005170	1	PRPT55 - Portique pied de lit Stryker SV1 - SV2	PRPT55 - Foot porch bar for Stryker SV1 - SV2
4655000071	1	B240 - Barre octogonale de 240cm	B240 - Plain bar of 240cm (94")
4655003398	2	DMEP - Double mâchoire pivotante	DMEP - Swivel double clamp
4655003067	1	T40E - Tube à mâchoire fixe de 40cm	T40E - Fixed clamp bar of 40cm (16")
4655000948	1	P753 - Étrier de traction	P753 - Handgrip
4655000617	7	PPO - Porte poulie ouverte	PPO - Open pulley with clamp
4635001304	3	SP01 - Sphère de protection	SP01 - Wall protection cushion

Construit à partir d'un matériau qui allie la résistance et la légèreté, le Système de Traction KINETEC® est d'un emploi facile pour tout le personnel soignant.

Kinetec SAS garantit son Système de Traction, 2 ans contre tous vices de fabrication, à partir de la date d'achat par le consommateur.

Dispositif médical de classe I conforme au règlement (UE) 2017/745.

Demandez notre catalogue complet d'accessoires et de kits de traction (poids, corde, attelle de positionnement, etc...)

Designed from a material which combines resistance and lightness, the KINETEC® TRACTION System is easy to use for all medical staff.

Kinetec SAS guarantees its Traction System 2 years against all defects of manufacture, starting from the date of purchase by the consumer.

Medical device of class I conforms to the regulation (EU) 2017/745.

Ask for our complete catalogue of accessories and traction kits (weights, cord, positioning splints, and more...)